



## Top 5 Care Tips for Tangible® Hydra-PEG® Coated Lenses

Coating your custom contact lenses with Tangible® Hydra-PEG® makes them feel wet, slippery and comfortable in the eye, and helps them stay clear and free of fogging.

Tangible Hydra-PEG coating is designed to improve comfort for the prescribed life of your lenses, but can be damaged if mishandled. Follow these tips to ensure that you protect your Tangible Hydra-PEG coated lenses.



### #1 Be careful your lenses don't slip off your fingers!

Lenses coated in Tangible Hydra-PEG will feel extra slippery, which helps the eyelids glide smoothly over your lenses. Take care when removing your coated lenses from their packaging or their contact lens case, and block your drain so you don't lose your lenses.



### #2 Use the correct cleaner.

**Approved cleaners for coated lenses are: Tangible® Clean, Boston Simplus, ClearCare, ClearCare HydraGlyde, Unique pH.** Always follow the instructions for use closely. If you have any questions, talk to your doctor.

Other contact lens cleaning solutions may damage the coating and shorten its benefits. **Do not use:** Progent, Boston Advanced Cleaner, Boston Laboratory Lens Cleaner, Boston One Step Liquid Enzymatic Cleaner, Optimum Extra Strength Cleaner, Optimum CDS, Quattro.

If you need to use artificial tears or rewetting drops, only use those approved for contact lenses.



### #3 Do not rinse lenses in tap water.

Do not rinse lenses in tap water. Only rinse your coated lenses with approved solutions or sterile saline. Do not use tap water, as it may wear away the coating.

Tap water also promotes the growth of bacteria and harmful pathogens on your contacts, which can cause eye infections. To protect against infection, rinse your contact lens case with approved solutions or sterile saline every day and allow the case to dry between uses.



### #4 Always store lenses in contact lens solution. Never store them dry.

Store your coated lenses with plenty of solution. If you leave them in your case dry, this can dry out and damage the special coating. If you need to store your lenses longer than overnight, make sure to change your contact lens solution monthly.



### #5 Last but not least, if you have any questions about caring for your lenses, talk to your doctor.

Your doctor knows your eye health history and is best equipped to recommend products, advise you about your contact lenses, and answer any questions.

**Enjoy your Tangible Hydra-PEG experience!**