



PROFESSIONAL SOLVENTS AND POLISH

- R1575** **Alox 721 Polishing Powder (1 lb. Tub)**
DESCRIPTION: Also called X-Pal. A white powder that needs to be mixed with water to make polish. Standard mixture is 12 parts water to one part 721 powder. Used for cleaning and polishing and all modifications on gas permeable or PMMA lenses.
- R1400** **Boston Laboratory Lens Cleaner (4 oz.)**
DESCRIPTION: A wax solvent intended to be used only in the practitioners office and not dispensed to patients. This solvent will remove deposits and waxy build-up on gas permeable lenses and should only be used once or twice on an individual lens. Damage may result if used too frequently or if a lens is allowed to soak in it for an extended period of time.
- R1416** **Boston Laboratory Lens Cleaner (16 oz.)**
DESCRIPTION: A wax solvent intended to be used only in the practitioners office and not dispensed to patients. This solvent will remove deposits and waxy build-up on gas permeable lenses and should only be used once or twice on an individual lens. Damage may result if used too frequently or if a lens is allowed to soak in it for an extended period of time.
- R1500** **Boston Professional Cleaning Polish (2 oz.)**
DESCRIPTION: Polishing compound for gas permeable contact lenses. Used for cleaning, polishing and all modifications on gas permeable or PMMA lenses. Pre-suspended; ready to use from the bottle.
- R1451** **Fluoro-Solve Wax Solvent (4 oz.)**
DESCRIPTION: A wax solvent intended to be used only in the practitioners office and not dispensed to patients. This solvent will remove deposits and waxy build-up, adhesive tape residue, body oils, and cosmetics on gas permeable lenses. When using solvents to clean lenses, it is best to minimize solvent-to-lens contact time (not to exceed 5 minutes).
- R1730** **Menicon Progent Starter *PROFESSIONAL IN-OFFICE USE ONLY!***
DESCRIPTION: Treatment for 1 pair of GP lenses plus reusable vial.
- R1732** **Menicon Progent 7-Pair Treatment Pack**
DESCRIPTION: Treatment for 7 pairs of GP lenses. Vial NOT included.
PROFESSIONAL IN-OFFICE USE ONLY!
- R1552** **Sil-02 Care Polish (32 oz.)**
DESCRIPTION: Polishing compound for gas permeable contact lenses. Used for cleaning, polishing and all modifications on gas permeable or PMMA lenses. Pre-suspended; ready to use from the bottle.
- R1550** **Sil-02 Care Polish (4 oz.)**
DESCRIPTION: Polishing compound for gas permeable contact lenses. Used for cleaning, polishing and all modifications on gas permeable or PMMA lenses. Pre-suspended; ready to use from the bottle.



R1730