

# Boston EO® ENVISION® Fitting Card



**Select  
9.6 mm  
Initial  
Diameter**

(Available in 9.3–10.3 mm)

**Determine  
Flat K  
(in Diopters)**

39.75	40.50	41.25	42.00	42.75	43.50	44.25	45.00	45.75	46.50	47.25	48.00
to	to	to	to	to	to	to	to	to	to	to	and
40.25	41.00	41.75	42.50	43.25	44.00	44.75	45.50	46.25	47.00	47.75	above



**Recommended  
Base Curve**

(Available in 8.3–7.0 mm)

**Corneal  
Astigmatism  
≤1.50D**

**Corneal  
Astigmatism  
>1.50D**

8.3	8.2	8.1	8.0	7.9	7.8	7.7	7.6	7.5	7.4	7.3
8.2	8.1	8.0	7.9	7.8	7.7	7.6	7.5	7.4	7.3	7.2

Trial Fitting Recommended



**Suggested  
Center  
Thickness**

**Lens Power**

**Thickness**

## Minus Lens Center Thickness

<b>Plano</b>	<b>-1.00</b>	<b>-2.00</b>	<b>-3.00</b>	<b>-4.00</b>	<b>-5.00</b>	<b>-6.00</b>	<b>-7.00</b>	<b>-8.00</b>
.18	.17	.16	.15	.14	.13	.13	.13	.13

Important: This lens has a clinically precise biaspheric back surface which already includes proper peripheral flattening to optimize the corneal fit.

**NO MODIFICATION SHOULD BE MADE TO THE BACK SURFACE, AS IT WOULD COMPROMISE CLINICAL PERFORMANCE.**