



BOSTON GAS PERMEABLE LENS CARE SYSTEMS

R1005

Boston Advance Comfort Formula Care System - #7470

Comes in quantities of one dozen.

INCLUDES: Boston Advance Cleaner (30 ml.)
 Boston Advance Comfort Formula Conditioning Solution (3.5 oz.)
 Boston One Step Liquid Enzyme
 Deep well lens case
 Patient instruction booklet

R1001

Boston Original Care System - #7475

Comes in quantities of one dozen.

INCLUDES: Boston Cleaner (30 ml.)
 Boston Conditioning Solution (3.5 oz.)
 Boston One Step Liquid Enzyme
 Deep well lens case
 Patient instruction booklet

R1050

Boston Ortho-K Care System - #7334

Comes in quantities of six.

INCLUDES: Boston Simplus Multi-Action Solution (3.5 oz.)
 Boston Rewetting Solution (10 ml.)
 DMV Lens Remover
 Deep well lens case
 Travel Bag

R1050



R1007

Boston Simplus Multi-Action Care System - #7333

Comes in quantities of one dozen.

INCLUDES: Boston Simplus Multi-Action Solution (3.5 oz.)
 Boston Rewetting Drops (5 ml.)
 Deep well lens case
 Patient instruction booklet

BOSTON GAS PERMEABLE LENS STARTER KITS

R1105

Boston Advance Comfort Formula Starter Kit - #7469

Comes in quantities of 24.

INCLUDES: Boston Advance Cleaner (10 ml.)
 Boston Advance Comfort Formula Conditioning Solution (1 oz.)
 Deep well lens case
 Patient instruction booklet

R1100

Boston Original Starter Kit - #7474

Comes in quantities of 24.

INCLUDES: Boston Cleaner (10 ml.)
 Boston Conditioning Solution (1 oz.)
 Deep well lens case
 Patient instruction booklet

R1107

Boston Simplus Multi-Action Starter Kit - #7332

Comes in quantities of 24.

INCLUDES: Boston Simplus Multi-Action Solution (1 oz.)
 Deep well lens case
 Patient instruction booklet



R1300

R1241



R1205

R1200

R1230

BOSTON GAS PERMEABLE LENS CLEANERS

- R1300 Boston Advance Cleaner (30 ml.)**
 Same solution and size that is included in the Boston Advance Comfort Formula Care System, and Boston Advance Comfort Formula Starter Kit (10ml bottle).
DESCRIPTION: Patient daily cleaner for gas permeable lenses.
NOT TO BE PLACED DIRECTLY IN THE EYE.
- R1241 Boston One Step Liquid Enzymatic Cleaner (4 week supply sample)**
 Comes in quantities of one dozen.
DESCRIPTION: Patient weekly cleaner for gas permeable lenses.
NOT TO BE PLACED DIRECTLY IN THE EYE.

BOSTON GAS PERMEABLE CONDITIONING/WETTING SOLUTIONS

- R1205 Boston Advance Comfort Formula Conditioning Solution (3.5 oz.)**
DESCRIPTION: Same solution that is included in the Boston Advance Comfort Formula Care System for the wetting and storage of gas permeable lenses. Preserved with PAPB, CHG and EDS. It is an effective agent that totally disinfects gas permeable lenses and kills all micro organisms when lenses are soaked for at least four hours.
- R1200 Boston Original Conditioning Solution (3.5 oz.)**
DESCRIPTION: Same solution that is included in the Boston Original Care Systems for wetting and storage of gas permeable lenses. Preserved with chlorhexidine. This preservative is not as effective as PAPB (polyaminopyl biguanide) but is a good alternative for those patients who react to PAPB.
- R1230 Boston Rewetting Drops (10 ml.)**
 Available in individual bottles or quantities of one dozen.
R1231 - 5 ml. samples also available in limited quantities - 20 samples/box.
DESCRIPTION: Rewetting solution to be used by patients during the wearing time to rewet lenses while they are still in the eye. Also preserved with PAPB (polyaminopropyl biguanide).