

LOANER SETS
All available for use up to 30 days

R1800

AKS Standard Cone Loaner Fitting Set

Diagnostic 13 lens keratoconic set with AKS design.

INCLUDES: Base Curves: 7.20 - 6.00 in .10mm steps
Diameter: 8.5

R1805

AKS Extended Cone Loaner Fitting Set

Diagnostic 13 lens keratoconic set with AKS design.

INCLUDES: Base Curves: 6.20 - 5.00 in .10mm steps
Diameter: 8.5

R1810

Boston Envision Loaner Fitting Set

Diagnostic 14 lens set.

INCLUDES: Base Curves: 8.30 - 7.00 in .10mm steps
Diameter: 9.6 (8.30 - 7.30 Base Curve), 9.3 (7.20 - 7.00 Base Curve)

R1815

Reverse Geometry Loaner Fitting Set

Diagnostic 11 lens reverse geometry set. Green.