

**TRIAL LENS SETS**  
*(Must Ship Traceable Method; All Sales Final)*

- R1609**            **Achievement Trial Lens Set (Boston ES\*)**  
Diagnostic 11 lens Achievement set with wet case.  
*INCLUDES:*    Base Curves: 7.30-8.30; Diameter: 9.40; Power: -3.00
- R1622**            **AKS Standard Diagnostic Fitting System**  
Diagnostic 13 lens set with storage case.  
*INCLUDES:*    Base Curves: 7.20 - 6.00; Diameter: 8.5; Powers: -6.00 to -12.00
- R1624**            **AKS Extended Range Diagnostic Fitting System**  
Diagnostic 13 lens set with storage case.  
*INCLUDES:*    Base Curves: 6.20 - 5.00; Diameter: 8.5; Powers: -11.00 to -17.00
- R1610**            **Artcon Trial Lens Set (Boston ES\*)**  
Diagnostic 11 lens Artcon set with wet case.  
*INCLUDES:*    Base Curves: 7.30-8.30; Diameter: 9.00 or 9.50; Power: -3.00
- R1600**            **Custom Trial Lens Set (any material)**  
Diagnostic 10 lens set with wet case; Power: -3.00
- R1601**            **Custom Trial Lens Set (any material)**  
Diagnostic 10 lens set without wet case; Power: -3.00
- R1606**            **Envision Trial Lens Set (Boston ES\*)**  
Diagnostic 14 lens set.  
*INCLUDES:*    Base Curves: 8.30 - 7.00 in .10mm steps  
                  Diameter: 9.6 (8.30 - 7.30 Base Curve), 9.3 (7.20 - 7.00 Base Curve)
- R1605MAG**        **MagniClearplus Trial Lens Set (Boston ES\*)**  
Diagnostic 10 lens MagniClearplus set with wet case.  
*INCLUDES:*    Base Curves: 7.30-8.20; Diameter: 9.50; Power: -3.00; Add Power: +2.00
- R1611**            **Boston Multivision Trial Lens Set (Boston ES)**  
Diagnostic 10 lens Boston Multivision set with wet case.
- R1615**            **Renovation Trial Lens Set (Boston ES)**  
Diagnostic 12 lens Renovation multifocal set.  
*INCLUDES:*    Base Curves: 8.30-7.20; Diameters: 9.50 (8.30-7.50 BC's) 9.20 (7.40-7.20 BC's)  
                  Power Dist: -3.00; Power Near: +2.50
- R1614**            **Thinsite 2 Trial Lens Set (Boston XO<sub>2</sub>)**  
Diagnostic 14 lens Thinsite 2 set with wet case.  
*INCLUDES:*    Base Curves: 7.30-7.90; Diameter: 9.00; Power: -3.00  
*OR:*            Base Curves: 7.70-8.30; Diameter: 9.50; Power: -3.00
- R1696**            **SoClear Standard Trial Lens Set (Boston XO)**  
Diagnostic 12 lens SoClear Standard set with wet case.  
*INCLUDES:*    Base Curves: 46.00-40.50 in .50D steps; Diameter: 14.00; Power: -3.00
- R1697**            **SoClear Keratoconic Trial Lens Set (Boston XO)**  
Diagnostic 12 lens SoClear Keratoconic set with wet case.  
*INCLUDES:*    Base Curves: 58.00-47.00 in 1.00D steps; Diameter: 14.00; Power: -3.00D at 47.00 to -14.00D at 58.00
- R1698**            **SoClear Individual Trial Lens**

(\*Boston ES is standard, may ship in Boston EO in backorder situations.)