

Item	Variables	KM Range	KM code	Power	Code	
Material	HDS100	40.00	A	04	-3.00	11.2 dia
Lens type	RG4	40.25	B	04	-3.00	
OAD (mm)	11.2 - 10.4	40.50	C	04	-3.00	
Edge Thickness	0.12	40.75	D	04	-3.00	
Index	1.442	41.00	E	04	-3.00	
CT	0.2	41.25	F	04	-3.00	
PWR	1.25	41.50	G	04	-3.00	
BOZ	6 mm	41.75	H	04	-3.00	
Laser	Marked	42.00	I	04	-3.00	
		42.25	J	04	-3.00	10.8 dia
		42.50	K	04	-3.00	
		42.75	L	04	-3.00	
		43.00	M	04	-3.00	
		43.25	N	04	-3.00	
		43.50	O	04	-3.00	10.4 dia
		43.75	P	04	-3.00	
		44.00	Q	04	-3.00	
		44.25	R	04	-3.00	
		44.50	S	04	-3.00	
		44.75	T	04	-3.00	
		45.00	U	04	-3.00	
		45.25	V	04	-3.00	
		45.50	W	04	-3.00	
		45.75	X	04	-3.00	
		46.00	Y	04	-3.00	
		46.25	Z	04	-3.00	
		46.50	ZA	04	-3.00	
		46.75	ZB	04	-3.00	
		47.00	ZC	04	-3.00	
		47.25	ZD	04	-3.00	
		47.50	ZE	04	-3.00	