



For Comfort and Clarity, the Choice is Clear — SoClear!

## SoClear<sup>®</sup> ASPHERICAL CONE

The comfort and stability achieved with the SoClear Aspheric Design gives the patient great visual acuity and comfort that surpasses that of standard RGP's for the vast majority of patients.

Variable aspheric back surface specially designed for moderate to advanced Keratoconus. Straight forward fitting methodology allows practitioners to achieve great results on hard to fit patients.

Provides Comfort and Stability for All Stages of:

Nipple Keratoconus

Oval Keratoconus

Fitting Set Base Curve range from 5.82 mm to 7.18 mm.

A Design from the SoClear<sup>®</sup> Family of Lenses  
***Clear, Comfortable, Healthy Vision***

### ART OPTICAL GUARANTEE:

Our worry-free approach to fitting GP lenses provides two refits per eye and full cancellation privileges for 120 days.



ORDERING 800.253.9364    CONSULTATION 800.566.8001    ONLINE [www.artoptical.com](http://www.artoptical.com)

Bausch & Lomb and Boston are registered trademarks of Bausch & Lomb Incorporated. 2009