



For Comfort and Clarity, the Choice is Clear — SoClear!

So₂Clear®
STANDARD

So₂Clear®
ASPHERICAL CONE

So₂Clear®
PROGRESSIVE

A lens design that is comfortable and stable, and can be used on a variety of indications. An essential tool that has broad application in the standard design range.

Stable Design Provides Unsurpassed Comfort and Crisp Vision for:

Excellent Astigmatism Correction

Difficult Single Vision Patients

Oval Keratoconus

Irregular Cornea

Pellucid Marginal Degeneration

Corneal Ectasia

Variable aspheric back surface specially designed for moderate to advanced Keratoconus. Straight forward fitting methodology allows practitioners to achieve great results on hard to fit patients.

Provides Comfort and Stability for All Stages of:

Nipple Keratoconus

Oval Keratoconus

The SoClear Progressive is a true revolution for today's Presbyopes. Imagine it — offering Presbyopes a multifocal lens that combines the best attributes of soft and GP multifocal lenses in one!

A True Revolution to Satisfy Presbyopes

Customizable Adds and Zone Sizes

Three designs from the SoClear® Family of Lenses
Clear, Comfortable, Healthy Vision

ART OPTICAL GUARANTEE:

Our worry-free approach to fitting GP lenses provides two refits per eye and full cancellation privileges for 120 days.

ARTOptical
contact lens, inc.

Bausch & Lomb
Boston.

ORDERING 800.253.9364 CONSULTATION 800.566.8001 ONLINE www.artoptical.com

Bausch & Lomb and Boston are registered trademarks of Bausch & Lomb Incorporated. 2009