

## TRIAL LENS SETS

- R1622**  
**R1800** (LOAN) **AKS Standard Cone Diagnostic Fitting System**  
Diagnostic 13 lens set with storage case.  
*INCLUDES:* Base Curves: 7.20 - 6.00 in .10mm steps  
Diameter: 8.5; Powers: -6.00 to -12.00
- R1624**  
**R1805** (LOAN) **AKS Extended Range Cone Diagnostic Fitting System**  
Diagnostic 13 lens set with storage case.  
*INCLUDES:* Base Curves: 6.20 - 5.00 in .10mm steps  
Diameter: 8.5; Powers: -11.00 to -17.00
- R1606**  
**R1810** (LOAN) **AKS Boston Envision Trial Lens Set (Boston ES)**  
(Used for Keratoconus, Transplants, Post RK, Post Lasik)  
Diagnostic 14 lens set.  
*INCLUDES:* Base Curves: 8.30 - 7.00 in .10mm steps  
Diameter: 9.6 (8.30 - 7.30 Base Curve), 9.3 (7.20 - 7.00 Base Curve)
- R1815** (LOAN) **Reverse Geometry Loaner Fitting Set in Boston EO material**  
Diagnostic 11 lens reverse geometry set.
- R1609** **Achievement Trial Lens Set (Boston ES)**  
Diagnostic 11 lens Achievement set with wet case.  
*INCLUDES:* Base Curves: 8.30-7.30; Diameter: 9.40; Power: -3.00
- R1610** **Artcon Trial Lens Set (Boston ES)**  
Diagnostic 11 lens Artcon set with wet case.  
*INCLUDES:* Base Curves: 8.30-7.30; Diameter: 9.00 or 9.50; Power: -3.00
- R1600** **Custom Trial Lens Set (any material)**  
Diagnostic 10 lens set with wet case; Power: -3.00
- R1601** **Custom Trial Lens Set (any material)**  
Diagnostic 10 lens set without wet case; Power: -3.00
- R1605MAG** **MagniClearplus Trial Lens Set (Boston ES)**  
Diagnostic 10 lens MagniClearplus set with wet case.  
*INCLUDES:* Base Curves: 8.20-7.30; Diameter: 9.50; Power: -3.00; Add Power: +2.00
- R1611** **Boston Multivision Trial Lens Set (Boston ES)**  
Diagnostic 10 lens Boston Multivision set with wet case.
- R1615** **Renovation Trial Lens Set (Boston ES)**  
Diagnostic 12 lens Renovation multifocal set.  
*INCLUDES:* Base Curves: 8.30-7.20; Diameters: 9.50 (8.30-7.50 BC's) 9.20 (7.40-7.20 BC's)  
Power Dist: -3.00; Power Near: +2.50
- R1614** **Thinsite 2 Trial Lens Set (Boston XO<sub>2</sub>)**  
Diagnostic 14 lens Thinsite 2 set with wet case.  
*INCLUDES:* Base Curves: 7.90-7.30; Diameter: 9.00; Power: -3.00  
*OR:* Base Curves: 8.30-7.70; Diameter: 9.50; Power: -3.00
- R1696** **SoClear Standard Trial Lens Set (Boston XO)**  
Diagnostic 12 lens SoClear Standard set with wet case.  
*INCLUDES:* Base Curves: 46.00-40.50 in .50D steps; Diameter: 14.00; Power: -3.00
- R1697** **SoClear Keratoconic Trial Lens Set (Boston XO)**  
Diagnostic 12 lens SoClear Keratoconic set with wet case.  
*INCLUDES:* Base Curves: 58.00-47.00 in 1.00D steps  
Diameter: 14.00; Power: -3.00D at 47.00 to -14.00D at 58.00
- R1698** **SoClear Individual Trial Lens**

Purchased Trial Sets are not returnable; Trial Sets available for loan must be returned within 30 days to receive credit.  
Trial Sets must ship traceable method.