

HANDLING TINT OPTIONS BY GP LENS MATERIAL

<u>MATERIAL NAME</u>	<u>COLORS AVAILABLE</u>
Alberta S	Blue (w/ UV)
Boston II.....	Blue (w/o UV)
Boston IV.....	Blue (w/o UV)
Boston EO.....	w/UV Blue, Electric Blue, Ice Blue, Green, Grey, Brown (without UV available in Blue & Green only)
Boston Equalens.....	Blue (w/ UV)
Boston Equalens 2.....	w/UV Blue, (without UV available in Yellow and Red only)
Boston ES.....	w/UV Blue, Ice Blue, Green, Grey, Brown, Clear (without UV available in Blue only)
Boston XO.....	w/UV Blue, Ice Blue, Green, Violet (without UV available in Yellow and Red only)
Boston X02.....	w/UV Blue, Ice Blue, Green, Violet (without UV available in Blue, Yellow & Red only)
Flosi.....	Blue, Blue #2 & #3, Clear, Green, Green #2 & #3, Grey, Grey #3, Brown, Violet, Yellow (w/o UV)
Fluorperm 30.....	Blue, Majestic Blue, Green, Grey, Clear (with UV available in Blue, Green & Grey only)
Fluorperm 60.....	Blue, Crystal Blue, Green, Clear, Brown (with UV available in Blue & Green only)
Fluorperm 92.....	Blue, Green, Clear (with UV available in Blue & Green only)
Fluorperm 151.....	Blue (available with or without UV)
Fluorex 300.....	Blue, Aqua Blue, Green, Grey, Clear, Rose Brown (not available w/ UV)
Fluorex 500.....	Blue, Aqua Blue, Green, Clear, Rose Brown (not available w/ UV)
Fluorex 700.....	Blue, Aqua Blue, Grey, Clear, Rose Brown(not available w/ UV)
Hydro2.....	Soft Blue, Cool Green, Ocean Blue (not available w/ UV)
Hybrid FS.....	Blue, Green, Grey, Clear (not available w/ UV)
Menicon Z.....	Blue (not available w/ UV)
ONSI 56.....	Blue, Green, Grey, Onzure Blue (with UV available in Blue only)
Optacryl 60.....	Blue w/o UV
OP-3.....	Blue, Dark Blue, Green, Dark Green, Grey, Clear, Brown w/o UV
Optimum Classic.....	Blue, Glacier Blue, Green, Grey w/ UV
Optimum Comfort.....	Blue, Glacier Blue, Green, Grey, Forest Green, Brown w/ UV
Optimum Extra.....	Blue, Glacier Blue, Green, Grey, Clear w/ UV
Optimum Extreme.....	Blue, Green w/ UV
Optimum HR 1.51.....	Blue (not available with UV)
Optimum HR 1.53.....	Blue (not available with UV)
Paragon HDS.....	w/o UV Blue, Green, Crystal Blue, Forest Green (with UV available in Blue, Green & Crystal Blue only)
Paragon HDS 100.....	Blue, Green (not available w/ UV)
Paragon HDS HI 1.54.....	Blue, Green (w/ UV only)
Paragon Thin.....	Blue, Green (with UV available in Blue only)
Paraperm 02.....	Blue, Electric Blue, Green, Green #2, Clear (w/o UV only)
Paraperm EW.....	Blue, Green, Clear (w/o UV only)
PMMA.....	Blue, Blue #2, Blue #3, Green, Green#2, Green#3, Grey, Grey #2, Brown, Brown#2, Clear, Violet, Red, Pink, Pink #2, Yellow/Green (w/o UV only)
SA18.....	Clear, Blue, Blue #2, Blue #3, Brown, Green, Green #2, Green #3, Grey, Grey #2, Violet, Red (w/o UV only)
SA32.....	Clear, Blue, Blue #2, Blue #3, Brown, Green, Green #2, Green #3, Grey, Grey #2, Violet (w/o UV only)
SGP 1.....	Blue, Green, Grey, Brown (w/o UV only)
SGP 2.....	Blue, Green, Grey, Clear, Brown (w/o UV only)
SGP 3.....	Blue, Green, Clear (w/o UV only)
Tyro 97.....	Blue, Green, Blue/Green, Grey, Clear (available with or without UV)
X Chrom.....	Red (w/o UV only)

HANDLING TINT AVAILABILITY BY LENS DESIGN:

Achievement ES.....	w/ UV Blue, Ice Blue, Green, Grey, Brown, Clear (without UV available in Blue only)
Achievement XO.....	w/ UV Blue, Ice Blue, Green, Violet
ARTCON.....	w/ UV Blue, Ice Blue, Green, Grey, Brown, Clear
Bi-Expert Bifocal.....	ES (Blue, Ice Blue, Green, Grey, Brown, Clear, without UV available in Blue only)
.....	EO (Blue, Electric Blue, Ice Blue, Green, Grey, Brown, without UV available in Blue & Green only)
.....	XO/XO2 (Blue, Ice Blue, Green, Violet)
Boston Multivision.....	Blue (w/ UV)
Lifestyle/Lifestyle Marquis.....	Blue (w/o UV)
OrthoTool.....	Yellow (HDS 100) (w/o UV)
Paragon CRT/Z-CRT/Sure Fit/Dual Axis/.....	Green (HDS 100 or HDS), Blue (HDS, HDS 100 or Menicon Z) (w/o UV)
RG4.....	Green or Blue (HDS 100) (w/o UV)
So ₂ Clear.....	w/UV Ice Blue (XO fenestrated, XO2 fenestrated or no-fenestrated, Optimum Extra no-fenestrated), Clear (Optimum Extra no-fenestrated, HDS 100 fenestrated or no-fenestrated), Crystal Blue (HDS 100 fenestrated or no-fenestrated)
Thinsite.....	ES (Blue, Ice Blue, Green, Grey, Brown, Clear, w/o UV available in Blue only)
.....	EO (Blue, Electric Blue, Ice Blue, Green, Grey, Brown, w/o UV available in Blue & Green only)
Thinsite 2.....	XO2 (Blue, Ice Blue, Green, Violet), ES (Blue, Ice Blue, Green, Grey, Brown, Clear, w/o UV available in Blue only)
.....	EO (Blue, Electric Blue, Ice Blue, Green, Grey, Brown, w/o UV available in Blue & Green only)
Envision.....	ES/EO/XO/XO2/EQ11 (all colors but red & yellow)