

TABLE 6 HIGH REFRACTIVE INDEX MATERIALS				
PRODUCT	MANUFACTURER	DK	SPECIFIC GRAVITY	COLORS
Optimum HR 1.51	Contamac US, Inc.	50	1.04	Blue
Optimum HR 1.53	Contamac US, Inc.	26	1.04	Blue
Paragon HDS HI 1.54	Paragon Vision Sciences	30	1.12	Blue, Green

Table 6 High Refractive Index Materials

These materials, with refractive indexes from 1.51 to 1.54, will assist in decreasing lens mass in single vision designs and can enhance the add power effect in multifocal designs.

*Quoted Dk values are measured using Gas to Gas method unless *, which indicates ISO/Fatt method. ¹FDA approved for extended wear. Many materials are available with or without UV - ask a customer service representative for UV information on your material of choice.*

MPC-310-T6