

Artcon® - crisp visual acuity with predictable fitting

ARTCON® combines all the best features of the **Polycon® II** lens design with the high performance of **Boston ES®** material. The result is a thin lens that offers crisp visual acuity, predictable fitting, excellent surface characteristics, durability and oxygen delivery. Combined with Art Optical's precision manufacturing, superior service and expedited delivery, **ARTCON®** is the perfect replacement lens for your **Polycon® II** patients.

PRODUCT FEATURES & BENEFITS

<i>Feature</i>	<i>Advantage</i>	<i>Benefit</i>
ARTCON® Lens Design	Polycon® II parameters, including the ultra-thin design	Comfortable, easy to fit, easy to upgrade current wearers
Manufactured in Boston ES® (Lowest silicon containing lens)	Added stiffness and strength, enhanced wettability, deposit resistance, crisp optics	Excellent comfort, performance, greater physiological response and durability
Design & Material Compatibility	Synergistic performance	Predictable fitting, performance and comfort
Wide Range of Base Curves & Powers	Fits lens to patient, not patient to lens	Higher fitting success



Art Optical and Boston® Guarantee: Unlimited exchanges, full cancellation privileges, and breakage/wettability coverage for a full six months provides a worry free approach to GP lens fitting.

ARTCON®

PARAMETER RANGE

Diameter: 9.0 and 9.5

Center Thickness: .08mm @ -5.00

Base Curves: 7.00mm to 8.85mm in .05mm steps

Powers: -20.00 to +20.00 in .25 diopter steps

Colors: blue, ice blue, green, gray, brown

ARTOptical
contact lens, inc.

Driven by your GP challenges
focused on your GP success

toll-free ordering **800.253.9364**

consultation direct **800.566.8001**

toll-free fax **800.648.2272**

online **www.artoptical.com**