



Ampleye Scleral Standard Diagnostic Set - Item# R1602

9 Ampleye Dx Lenses with 16.5mm diameter and toric haptic. *INCLUDES:* Sagittal depths of 3800, 4000, 4200, 4400, 4600, 4800, 5000, 5200, & 5600 in storage case with fitting instructions and accessories - large DMV, small DMV, hand-held cobalt blue light, wrattan filter, & starter supply of non-preserved saline.

Ampleye Scleral 15.0/15.5 Standard Diagnostic Set - Item# R1603

4 Ampleye Dx Lenses with 15.0mm diameter and 5 Ampleye Dx Lenses with 15.5mm diameter and toric haptic. *INCLUDES:* Sagittal depths of 3400, 3800, 4200, 4600 in both diameters and 4800 in 15.5mm in storage case with fitting instructions and accessories - large DMV, small DMV, hand-held cobalt blue light, wrattan filter, and a starter supply of non-preserved saline.



Ampleye Scleral 15.0/15.5/16.5 Combo Diagnostic Set - Item# R1605

18 Ampleye Dx Lenses with 15.0/15.5/16.5mm diameter. *INCLUDES:* 1 lens of 3400, 3800, 4200, 4600 in both 15.0 & 15.5 diameters and 4800 in 15.5mm, as well as, one each of sag depths 3800, 4000, 4200, 4400, 4600, 4800, 5000, 5200, & 5600 in 16.5mm. Shipped in storage case with fitting instructions and accessories - large DMV, small DMV, hand-held cobalt blue light, wrattan filter, and starter supply of non-preserved saline.

Ampleye Scleral Diagnostic Set Extensions - Item# R1602EXT or R1603EXT

Sold as add-ons to Ampleye R1602 or R1603 9-lens sets. **R1602EXT** includes duplicate 16.5mm lenses with 3800, 4000, 4200, 4400, 4800 as well as a 5400sag, all with toric haptic, and spherical lenses of 4000, 4200 & 4400sag. **R1603EXT** includes (4) 15.0mm spherical lenses & (5) 15.5mm w/toric haptic in 3400, 3800, 4200, 4600sag, plus 4800sag in 15.5mm. Lenses will ship with an 18-lens foam insert that fits the original storage case.



Ampleye Individual Trial Lens - Item# R1602IRL - Available to order by diameter & sag for existing set-holders.



Ampleye Patient Starter Kit Item# R1130

INCLUDES: deep well lens case, large DMV for lens insertion, traditional DMV for lens removal, pocket mirror, starter supply of non-preserved saline, & instructions in clear carrying case.

Nutrifill - Item# R1128

Preservative free rinsing & insertion solution for scleral lenses. Formulated with 5 electrolytes to mimic natural tears. Limited sample quantities available to practitioners for S&H fee. Wholesale and patient supplies available @ nutrifill.com. Made in USA.



Menicon LacriPure Saline Solution Item# R1127

98 Count Box - 5ml vials. FDA approved non-preserved saline solution indicated for rinsing soft, hybrid, or GP lenses and for scleral lens insertion.

DMV Traditional RGP Lens Removers Item# R3100

Small suction cup for removing gas permeable lenses from the eye. Instructions included.



Progent by Menicon

For biweekly cleaning of scleral and corneal GP lenses. Cleans, loosens and removes surface deposits after a short soak without manual rubbing. *INCLUDES:* Progent & lens case.

DMV Scleral Cup - Item# R3130

Large diameter GP lens insertion tool.



EZI Scleral Lens Applicator Item# R3135

Scleral lens insertion ring allows for one finger application, lens self-positioning, and less air entrapment.



Item# R1730 - Single Treatment

Item# R1732 - 7 Treatment Pack

Progent by Menicon Scleral Lens Case Item# R1746

For Progent treatment of large diameter lenses (11-23mm).

