

Intelliwave® Conversion Guide for Coopervision® Hydrasoft® Patients

Lens Design CooperVision Hydrasoft Sphere/Toric	Material Methafilcon B 55%	Dk 18.80	Base Curve/Diameter 8.3, 8.6/14.2 8.3, 8.6, 8.9, 9.2/15.0	Powers +30 to -30D	Cylinder -0.50 to -10.00D	Axis 1°	Replacement Quarterly
Comparable Option Intelliwave Asphere/Toric	Material Acofilcon B 49%	Dk 15.89	Base Curve/Diameter 7.8 to 9.5/13.0 to 15.5 Order Exact Duplicate	Powers* +30 to -30D	Cylinder -0.50 to -12.00D	Axis 1° to 180°	Replacement Semi-Annual or Quarterly
Higher Dk Option Intelliwave Asphere/Toric	Material Eprofilcon A 74%	Dk 60	Base Curve/Diameter 7.8 to 9.5/13.0 to 15.5 Steepen Base Curve by .20mm	Powers* +30 to -30D	Cylinder -0.50 to -12.00D	Axis 1° to 180°	Replacement Quarterly

Need a Multifocal?

Intelliwave Multifocal & Multifocal Toric lenses are available in Center Near or Center Distance design options; Add Powers from +0.50 to +4.00D. Contact our professional consultation team at **800.566.8001** with additional questions - they are always happy to assist!

*powers greater than +/-20.00 are considered off-label ©AOCLI 2022 Intelliwave is a registered trademark of Art Optical Contact Lens, Inc. - Any other product/brand names are trademarks of the respective owners.

ARTOptical
contact lens, inc.

proudly featuring quality lens materials by

Contamac®