

SoClear[®]

Corneal-Scleral Contact Lens Family

Three designs from the SoClear[®] Family of Lenses

Standard

A lens design that is comfortable and stable, and can be used on a variety of indications. An essential tool that has broad application in the standard design range.

Stable Design Provides Unsurpassed Comfort and Crisp Vision for:

- Excellent Astigmatism Correction
- Difficult Single Vision Patients
- Oval Keratoconus
- Irregular Cornea
- Pellucid Marginal Degeneration
- Corneal Ectasia

Fitting Set Base Curve Range from 7.34 mm to 8.33 mm

Aspherical Cone

Variable aspheric back surface specially designed for moderate to advanced Keratoconus. Straight forward fitting methodology allows practitioners to achieve great results on hard to fit patients.

Provides Comfort and Stability for All Stages of:

- Nipple Keratoconus
- Oval Keratoconus

Fitting Set Base Curve Range from 5.82 mm to 7.18 mm

Progressive

The SoClear Progressive is a true revolution for today's Presbyopes. Imagine it - offering Presbyopes a multifocal lens that combines the best attributes of soft and GP multifocal lenses in one!

- A True Revolution to Satisfy Presbyopes
- Customizable Adds and Zone Sizes

Fitting Set Base Curve Range from 7.34 mm to 8.33 mm

For Comfort and Clarity,
the choice is Clear - SoClear!
Clear, Comfortable, Healthy Vision

One Warranty, No Worries

Unlimited Exchanges! No Lens Returns Required*! No Material Exchange Fees! Our worry-free fitting warranty has you covered for 120 days from the initial order date.

ARTOptical
contact lens, inc.

toll-free ordering **800.253.9364**
consultation direct **800.566.8001**
toll-free fax **800.648.2272**
online **www.artoptical.com**

*Exceptions apply for cancellation and non-parameter or Rx related exchanges. Policy subject to change.