

SOCLEAR TRIAL LENS SETS

- R1696** **SoClear Standard Trial Lens Set (Boston XO)**
Diagnostic 12 lens SoClear Standard set with wet case.
INCLUDES: Base Curves: 46.00-40.50 in .50D steps; Diameter: 14.00; Power: -3.00
- R1697** **SoClear Keratoconic Trial Lens Set (Boston XO)**
Diagnostic 12 lens SoClear Keratoconic set with wet case.
INCLUDES: Base Curves: 58.00-47.00 in 1.00D steps
Diameter: 14.00; Power: -3.00D at 47.00 to -14.00D at 58.00
- R1698** **SoClear Individual Trial Lens**
- R1696F** **SoClear Standard Fenestrated Trial Set (Boston XO)**
Diagnostic 12 lens SoClear Standard set with wet case.
INCLUDES: Base Curves: 46.00-40.50 in .50D steps; Diameter: 14.00; Power: -3.00
- R1697F** **SoClear Keratoconic Fenestrated Trial Set (Boston XO)**
Diagnostic 12 lens SoClear Keratoconic set with wet case.
INCLUDES: Base Curves: 58.00-47.00 in 1.00D steps
Diameter: 14.00; Power: -3.00D at 47.00 to -14.00D at 58.00
- R1698F** **SoClear Individual Fenestrated Trial Lens**

Purchased Trial Sets are not returnable. Trial Sets must ship traceable method.